

## User Tips & Reminders for the *Celluma*:

- Do not lean into or place full body weight on the device
- Please reference the user manual for protocols, contraindications and maintenance of the device
- The device should be stored flat when travelling or when not in use
- Place as close as possible to area you are treating
- For facial treatments, remove make-up as many formulations contain minerals which deflect the light – same applies with sunblock
- Place device close to clean dry skin (move clothing away)
- It is OK to use *Celluma* daily
- Be creative with placement – there is usually more than one way to place the device over any given body part (but always keep your client comfortable and unstrained)
- Place legs on an ottoman
- Except for facial treatments it is often possible to seat clients so they can enjoy a *Celluma* treatment while reading or using a laptop/smart phone/iPad
- Return your warranty card: serial number can be found on the back of the white controller starting with a “C”
- Don’t forget *Celluma* for your pets. Animals seem to know what’s inherently good for them and respond really well to *Celluma*.

## How to Use *Celluma* Over the Face:

- To ensure ideal positioning, pre-form *Celluma* into a “C” shape.
- Alternatively, wrap around head – may be a bit uncomfortable though
- Regular tent – OK, but a bit too far from face
- Raising head on pillow or neck roll to move closer to diodes, then “tuck” ends under pillow
- Vertically – covering face, neck and décolleté

NOTE: If you purchased a *Celluma LITE* and wish to use it over the face you will need a pair of *Celluma LITE Facial Rests*.